E

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GB

Affections of the pelvis ring (e.g. slackness of the connective tissue of several pelvic bones) and in connection with it instability are a frequently occurring after-effect of the reorientation process of the hormonal household of people passing the age of 50. Females - not at east caused by pregnancy - are certainly more frequently affected than males.

Instructions for use:

Position the brace in the middle on the pelvis so that the upper edge of the brace is not higher than the upper iliac crest or protruding into the waist. Insert a finger in the finger loop provided on the fastener of the abdominal plate. Now pull the fastener evenly forwards.

Press the left half of the fastener onto the abdomen and place the right half of the fastener over the left half until you can attach it on the hook and loop pad on the left fastener half. When securing, withdraw the finger first from the left and then from the right finger loop.

Now grasp both tension straps in both hands and pull them simultaneously and steadily forwards until the desired force effect (stabilisation) is produced. Then secure both tension straps with the hook and loop fastener on the abdominal plate.

Indications: Symphysis insufficiency, symphysis rupture, structural loosening in the area of the iliosacral joints / pubic symphysis, after birth.

Contraindications: Peripheral arterial occlusive disease (PAOD), lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid, sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body.

Composition: 44 % polyamide, 25 % polyester, 25 % viscose, 4 % spandex, 2 % cotton

Important instructions:

No reutilisation – this medical device is meant for the care of one patient only.

In case of doctor's order to wear support during night, avoid hindering blood circulation. In case of numbness loosen support and remove if necessary. If complaints continue, please contact the doctor.

The product contains latex and can trigger allergic reactions.

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