

RANGE OF MOTION WALKER

Range of Motion walkers offer the ability to allow ankle motion or fix the ankle at a desired angle.

The range of motion joint adjusts in 10° increments, up to 40° in both dorsi and plantar flexion, and can either be completely locked out at any desired interval, or left with desired range and requires no additional hardware to set ROM.

Designed with wider rocker bottom soles to yield greater stability, the ROM walker is indicated for aiding acute trauma, moderate to severe ankle sprains and fractures and is ideal for patients who have had achilles tendon repairs.

ORDERING INFORMATION

PART #	SIZE	MENS SHOE	WOMENS SHOE
ROM1200	XSmall	up to 4	up to 5.5
ROM1400	Small	4.5-7	6-8
ROM1600	Medium	7.5-10.5	8.5-11.5
ROM1800	Large	10.5-12.5	11.5-13.5
ROM2000	XLarge	12.5 & up	13.5 & up



FOLLOW ÖSSUR ON



WWW.OSSUR.COM.AU



Össur Australia - OA & Injury Solutions
TEL 1300 123 268 TEL +61 3 8761 6408
FAX +61 3 9761 6067
infomelbourne@ossur.com

© ÖSSUR, 11. 2015 | P-800035-B